











### Instructions of use Ritex MEDICAL PROBE COVERS for sonography probes

Ritex protective covers for ultrasound probes undergo several quality tests and are characterised by very good skin compatibility. Each protective cover is individually electronically tested for leaks.

- Colour: transparent
- Surface structure: smooth
- Shape: cylindrical
- Material: natural rubber latex
- Coating: dry, powdered
- Condition: non-sterile
- Available content: 200 covers
- Available sizes: Ø 34 mm x 200 mm, Ø 28 mm x 200 mm

	Medical Device.
	Made of natural rubber latex which can provoke allergic reactions, including anaphylactic shock, in people who are allergic to latex and must not be used from those affected.
	Keep away from direct sunlight and heat.
	Keep in a dry place.
	Use only once as reuse poses a greater risk of damage or infection.
	Batch code.
	Expiry date.
	Manufacturer.
	Do not use any lubricants containing oil or grease (e.g. Vaseline, body lotion, etc.). Observe the product description of the ultrasound gel.
	Disposal: Dispose of with household waste.
	Reporting information: Report serious incidents involving the use of the product to the competent authority (location of user/patient).
CE	Tested using EN ISO 4074 as reference

Dated 23/05/2022